

## PCR using KAPA

- Dilute DNA template:
  1. 1 ng/ $\mu$ l if template is a PCR fragment
  2. 10 ng/ $\mu$ l if template is a plasmid
  3. 100 ng/ $\mu$ l if template is either a bacmid or genomic DNA
- Dilute primers to a final concentration in the reaction of 0.3  $\mu$ M:
  1. Take 1.5  $\mu$ l from each primer and add 12  $\mu$ l water – this gets you 10  $\mu$ M
  2. From the 10  $\mu$ M tube take 1.5  $\mu$ l into 50  $\mu$ l reaction
- Set up the reaction (for 50  $\mu$ l reaction):
  1. 1  $\mu$ l DNA template
  2. 1.5  $\mu$ l primer mix
  3. 22.5  $\mu$ l water
  4. 25  $\mu$ l 2X Kapa HiFi HotStart ReadyMix
- Thermocycler:
  1. 3 minutes at 95 °C
  2. 10-15 seconds at 98 °C
  3. 15-20 seconds at 60 °C (This is for annealing. Don't go below 60 °C . You can increase the temperature or use a gradient if you get non-specific amplification)
  4. 1 minute/kb template at 72 °C
  5. Repeat steps 2-4 for 29-34 cycles
  6. 5-10 minutes at 72 °C
  7. Keep at 12 °C