## **Bacmid transfection (P0)**

- 1. Count insect cells and plate 1 million cells per well in a 6-well plate. If the cells are too dense, dilute them to 0.5-1 cells/ml so that there is enough liquid to cover the bottom of the well.
- 2. Let incubate for 20 minutes in a humidity chamber with a wet paper towel
- 3. For each transfection mix:
  - a. Tube A 2 μgr bacmid + 100 μl grace media
  - b. Tube B 6 µl FuGene HD + 100 µl grace media
- 4. Let tubes incubate for 5 minutes
- 5. Mix tube A + tube B and incubate for 15 minutes
- 6. Remove media from cells
- 7. Add 800  $\mu$ l grace media to the cells, gently, on the wall of the well, so not to disturb the cells.
- 8. Add the transfection mix gently to the cells from above, drop-wise
- 9. Incubate for 4-5 hours
- 10. Add 2 ml of full media (ESF 921)
- 11. Incubate for 4-5 days. If using DH10EMBacVSV, monitor infection by fluorescence in the red channel
- 12. Harvest virus by spinning cells in a conical tube at 4000 rpm. Transfer supernatant to a new, sterile tube. Use all the P0 to infect 30 ml of insect cells at 1E6 cells/ml to generate P1. Incubate for 3-4 days
- 13. Make P2: use 3 different concentrations (typically 5  $\mu$ l, 50  $\mu$ l, and 500  $\mu$ l) of P1 to infect 200 ml of insect cells at 1E6 cells/ml. Incubate for 3 days.
- 14. Spin down at 4000 g for 10 minutes. Add 5% FBS and filter sterilize. This is your P2.